



# monoline

## Scintillation Detector with Integrally Mounted PMT

### ■ advantages

- Compact assembly
- Direct PMT to crystal mounting
- PMT is matched and tested with scintillator
- Consistent, superior energy resolution
- Resolutions for individual detectors can be guaranteed



In the Bicon *Monoline* assemblies, the photomultiplier tube is optically coupled directly to the scintillator. The scintillator is mounted in a container (usually of aluminum), and a mu-metal magnetic/light shield is fitted over the PMT. The scintillator container and mu-metal shield are sealed together to form a low-mass and light-tight housing for the detector.

This design usually yields better and more consistent pulse height resolution than others. Monolines are, therefore, the detectors of choice for spectroscopy and radioisotope assay.

### ■ options

- Scintillator containers of low-background stainless steel or copper
- Thin aluminum or beryllium radiation entrance windows
- PMTs selected for low background, premium resolution, fixed HV use, or gain matching
- Special flanges, mounting fixtures or other modifications
- Integrated, low-background voltage divider/preamp bases
- Square, hexagonal or other cross-sections
- Premium resolution assemblies

### ■ other configurations

- Waterproofed assemblies
- Ruggedized and high-temperature detectors
- Assemblies using thin crystals for low-energy gamma and x-ray detection

### ■ design notes

- Some configurations, especially those employing thin scintillators, may require light pipes between the scintillator and PMT. The light pipes diffuse the propagated light to overcome non-uniformity effects in the PMT photocathode.
- The detector package is hermetically sealed when NaI(Tl) or another hygroscopic scintillator is used.
- The maximum scintillator size is 127mm (5") in diameter or on the diagonal.




**BICRON**<sup>®</sup>

## popular configurations

### solid

Commonly used for gamma ray spectroscopy, radon canister counting, thyroid uptake measurements, health physics.

Model 3M3/3 shown




Typical Solid Crystal Monoline Models

Model	Crystal Size	PMT Size	Model	Crystal Size	PMT Size
1M1/1.5	1" x 1"	1.5"	2M2/2	2" x 2"	2"
1M2/2	1" x 2"	2"	3M2/3	3" x 2"	3"
1.5M.5/2	1.5" x .5"	2"	3M3/3	3" x 3"	3"
1.5M1/2	1.5" x 1"	2"	4M4/5	4" x 4"	5"
1.75M2/2	1.75" x 2"	2"	5M5/5	5" x 5"	5"

### end well

The most efficient configuration; for radioisotope assay, wipes, sample counting.

Model 2MW2/2 shown



Typical End Well Crystal Monoline Models


Model	Crystal Size	Well Size	PMT Size
2MW2/2	2" x 2"	see table	2"
3MW3/3	3" x 3"	see table	3"
3MW4/3	3" x 4"	1.25" x 3"	3"

Well Sizes	
2MW2/2	3MW3/3
.625" x 1.437"	.625" x 1.625"
.656" x 1.546"	.656" x 2.063"
.787" x 1.378"	1" x 2"
.75" x 1.437"	1.125" x 2"
1" x 1.546"	1.187" x 2.062"
1.125" x 1.437"	1.333" x 2"

### through-side well

An ideal configuration when space is limited; the second most efficient configuration; used in radioisotope assay and fuel rod monitoring.

Model 2MSW2/2 shown  
(.75" ID well)



Manufacturer reserves the right to alter specifications.



3005(01-97)

12345 Kinsman Road • Newbury, Ohio • 44065 • USA • Phone: (440) 564-2251 • Fax: (440) 564-8047 • <http://www.bicron.com>

SAINT-GOBAIN  
INDUSTRIAL CERAMICS

Saint-Gobain Industrial Ceramics, Inc.